

Install Courier and Select Lab

Courier is not compatible with Macs unless Bootcamp is installed.

1. Visit www.360courier.net

2. Click "Download"



3. Complete each field and click "Register & Download"

NAME	<input type="text" value="Joe Smith"/>
EMAIL	<input type="text" value="joe.smith@yahoo.com"/>
PRACTICE	<input type="text" value="Dr. Joe Smith, DDS PC"/>
PHONE	<input type="text" value="555-555-5555"/>
ADDRESS	<input type="text" value="123 Smiles Blvd. Dentaltown, AL 12345"/>

4. 360courier installer will download. Select "Run" and "Next" through the prompts.



5. An automatically generated email will be sent to the email address you provided during step 3. This email will include the Serial and Activation Code.

Dear Joe,

We are delighted for your interest in our valuable application 360Courier.

360Courier application is there to make it easy for you to share your DICOM images, 360dps plans and other attached files in a HIPPA compliant way, with your lab and also with your colleagues, submit, track and manage your implant planning and surgical guide orders.

360Courier is a valuable addition that complements the wonderful planning experience of your 360dps application.

Please download the latest version of courier from <http://www.360courier.net>

After installation, use the following serial and activation codes to start enjoying your experience.

Serial: `cd9843e782796a6c3ab4322791e0e082885`

Activation: `a897123340eae4ab095e09b77121ae4798`

Please be noted that your service contract is valid till: 2016/08/19

In case of lost or non working serials or any inquiries please inform 360Courier support team on: support@360courier.net

also, it is important to make sure that your personal information is correct as recorded in our database.

Your personal information are:

Name: **Joe Smith**

Email: joe.smith@yahoo.com

Mobile: 555-555-5555

for any mistaken data please contact us on: support@360courier.net


We hope you enjoy your 360Courier experience.

6. Copy and paste the codes into the 360courier field and click "Activate."




Selecting a Lab:


1. Select the name on the account. This is found in the top, right corner.




2. Select "Preferences"

A screenshot of a "User Preferences" dialog box. The dialog has three tabs: "Basic Information", "Practice Information", and "Preferences", with the "Preferences" tab selected and highlighted in red. The "Name" section contains two text input fields with the values "Joe" and "Smith". The "Contact Info" section contains two text input fields with the values "555-555-5555" and "joe.smith@yahoo.com". At the bottom right, there are "Cancel" and "Save Changes" buttons.


3. Click “360 Imaging” to choose a lab.



4. Select Oral Arts Dental Lab and click “OK”



5. Next, select “Practice Information”, complete the entire form, and click “Save changes.”



The image shows a 'User Preferences' dialog box with three tabs: 'Basic Information', 'Practice Information', and 'Preferences'. The 'Practice Information' tab is active. The form contains the following fields:

- Name: Dr. Joe Smith, DDS PC
- Profession: General Dentist (dropdown menu)
- Address section:
 - Address Line 1: 123 Smiles Blvd.
 - Address Line 2: (empty)
 - Country: USA (dropdown menu)
 - State: Alabama (dropdown menu)
 - City: Dentaltown
 - Zip: 12345

At the bottom right, there are two buttons: 'Cancel' and 'Save Changes'.

Now your lab of choice is selected and all uploaded files will go directly to that lab.