



Digital Impression with Scan Bodies Workflow

Modeless Screw-Retained Crown from Digital Scan:

Scan body required- Glidewell IOS kit from GL

Ti/Zr Abutment from Digital Scan:

Scan body required - DentaSwiss kit from Oral Arts

Protocol: *

- 1) Use appropriate scan body in place of your prep and log in to your digital scanning system accordingly
- 2) Please add notes on the digital RX notating the scan body kit used, which implant platform was used, etc.

*Can be used with all scanning systems, except Sirona. Oral Arts does not have a digital implant workflow for Sirona.

Please call Oral Arts directly with any questions:

Anthony "AC" Smithwick, Digital Crown and Bridge Supervisor:

1-800- 354-2075 Ext. 7122

Josh Cook, CDT, Implant Department Supervisor:

1-800- 354-2075 Ext. 7145

Implant Department:

1-800- 354-2075 Ext. 7178