

Model-Less v.s. Printed Model Crown Options

✓	Restoration type	Bite Scan	Posterior/Anterior	Number of Units
Model-Less Cases Accepted	Monolithic crowns (E.max, Element-Z, or Ultra-Z) and CAD/CAM Oral Temps	Accurate bite scan required	Posterior preferred, anterior accepted	Up to 2 consecutive single units, 3 unit bridges
Printed Models Recommended			Anterior restorations, enhanced esthetics achieved with model	
Printed Model Required	Multi-layered restorations, veneers, onlays, PFM, zirconia lingual, Forte, and Select	If drastic changes are required to the bite, we will need to print models and hand articulate		3 or more consecutive single units, 4 unit bridge or larger