

# Model-Less v.s. Printed Model Crown Options

✓	Restoration type	Bite Scan	Posterior/Anterior	Number of Units
<b>Model-Less Cases Accepted</b>	Monolithic crowns (E.max, Element-Z, or Ultra-Z) and CAD/CAM Oral Temps	Accurate bite scan required	Posterior preferred, anterior accepted	Up to 2 consecutive single units, 3 unit bridges
<b>Printed Models Recommended</b>			Anterior restorations, enhanced esthetics achieved with model	
<b>Printed Model Required</b>	Multi-layered restorations, veneers, onlays, PFM, zirconia lingual, Forte, gold crowns, and Select	If drastic changes are required to the bite, we will need to print models and hand articulate		3 or more consecutive single units, 4 unit bridge or larger